



Calibration Services to **GPX1500** Vial

A calibration package allowing an efficient calibration of the GPX1500 Vial instrument, that can be done remotely without having to send the instrument back and forth.

A Reference Standard set is a prerequisite to perform the Calibration Service. A Reference Standard are stoppered vials or flame sealed vials or ampoules certified by Gasporox with known concentrations.

With the Reference Standards the customer generate measurement with its instrument with the support of a manual.

The measured data is sent to Gasporox and verified by Gasporox support. After analyzing the data, a calibration certificate for the instrument with a validity of 12 months is issued and sent to the customer.

Supported Instruments

- GPX1500 vial base unit
- Format parts

Benefits

- No need to send in the instrument
- Instrument availability
- Fast calibration procedure
- Easy process – only three concentrations needed